MEDICAL PRODUCT CATALOGUE

2016



## **CONTENTS**

Diabetes

In - Patient Care

Blood Pressure Monitoring System

**Fabric Orthopedics** 

**Metal Orthopedics** 


**Foot Health** 

Orthopedic Equipments

**Respiratory Care** 

**Personal Care Products** 

Disposables





## DIABET



### **FORA**



#### **BLOOD GLUCOSE MONITORING SYSTEM** FORA G20

- No code
- Strip Ejection
- Function
- Measurement Range: 20 ~ 600 mg/dL
- (1.1~33.3mmol/L)
- Daily alarm: 4
- Test Strip: FORA G20 test strips
- SampleVolume: 0.7 µL
- Reaction Time: 7 seconds
- Coefficient of Variation: Less than 5%
- **GMemory Capacity:** 450 results with date & time
- Dimensions: 87 x 51 x 14 mm
- Weight: 49 g (batteries excluded)
- Power Source: One CR2032 Battery

# **FORA**



#### **BLOOD GLUCOSE STRIP - FORA G20**

- No codina
- Blood test volume (uL): 0.7
- Reaction time (seconds): 7
- Measurement range: 20 ~ 600 mg/dL (1.1~33.3mmol/L)
- Enzyme type: god
- Hematotocrit range: 20 to 60 %
- Memory set: 450
- Strip ejection: Yes
- Daily alarm: 4
- Alternative site testing: Yes
- Ketone warning:
- **Day average:** 7/14/21/28/60/90 days
- Data output: RS -232
- Power saving (seconds): 180
- Power source: 1 x CR 2032
- Dimension (Inch): 3.42" x 2.00" x 0.55"
- Weight (lb) without batteries: 0.11
- Unit per box: 50 strips





### **e**Bsensor



#### **BLOOD GLUCOSE MONITORING SYSTEM** eBsensor

- Measuring Unit: mg/dL or mmol/L
- **Size:** 87 x 60 x 21 mm
- Memory Capacity: 180 blood glucose results
- System Operating Temperature:  $10{\sim}40\,^{\circ}\text{C}$
- System Operating Humidity: < 85%
- Meter Storage Temperature: 0~50 °C
- Meter Storage Humidity: < 95%
- Power Supply: Two 1.5-Volt (AAA) batteries
- Weight: 75g







### **e**Bsensor



## **eB**fine



### **FORAFINE**



### **e**Bsensor



#### **BLOOD GLUCOSE STRIP - eBsensor**

- Blood Sample: Capillary Whole Blood
- Blood Volume: 2.5 µl • Hct Range: 20 % ~ 60 %
- Measuring Range: 20~600mg/dL (1.1

~33.3mmol/L)

• Measuring Time: 10 secs • Unit per box: 50 strips

#### **INSULIN PEN NEEDLES – eBfine**

- Triple facet bevel needle for a gentle and Min. puncture for Max. injection comfort
- Silicone lubricant coating on insulin pen for smooth
- Non-Toxic, non-Pyrogenic needle for guaranteed safety
- Optional size and fit major brands injection pens
- Unit per box: 100 unit

#### **INSULINE PEN NEEDLES – FORAFINE**

- Fits all major brands of injection pens
- Triple sharpened needles for easier penetration
- Silicone coated needle tips and thin wall technology for maximum injection comfort and safety
- · Easier and guicker insuline flow
- Unit per box: 100 unit

#### LANCET (100 pcs)

- Safe: Sterilized by gamma irradiation
- Tri-level ground cannula keeps good blood flow
- Universal: Optional size and fit for almost all lancet devices
- Single use only









ISO 13485 | ISO 9001



# **IN – PATIENT CARE**





## **IN - PATIENT CARE**



## Wicromed.



## Wicromed.



## Wicromed.



#### **ADULT DIAPER (MEDIUM) (8 PCS)**

- Unit per box: 8 units
- Box per package: 12 boxes
- Product dimensions: 70x84 cm
- Waist Sizes: 85x125 cm
- User Weight: 45-75 kg

#### **ADULT DIAPER (LARGE) (7 PCS)**

- Unit per box: 7 units
- Box per package: 12 boxes
- Product dimensions: 80x99 cm
- Waist Sizes: 100x150 cm
- User Weight: 75 kg and more

## ADULT DIAPER (MEDIUM) (30 PCS) (ECONOMIC)

- Unit per box: 30 units
- Box per package: 4 boxes
- Waist Sizes: 85x125 cm
- User Weight: 45-75 kg









## Wicromed.



## Wicromed.



## Wicromed.



#### **ADULT DIAPER (LARGE) (30 PCS)** (ECONOMIC)

• Unit per box: 30 units • Box per package: 4 boxes Waist Sizes: 100x150 cm • User Weight: 75 kg and more

- Textile backsheet With barriers
- Super Absorbent paritcles
- Extra Dryness Layers
- Suitable for sensitive skin
- Special type and frontal tape system
- Wetness Indicator
- Unit per box: 7 units • Waist Sizes: 100x150 cm • User Weight: 70 kg and more

### **ADULT DIAPER (SMALL) (30 PCS)**

- Textile backsheet
- With barriers
- Super Absorbent paritcles
- Extra Dryness Layers
- Suitable for sensitive skin
- Special type and frontal tape system
- Wetness Indicator
- Unit per box: 30 units • Waist Sizes: 50x85 cm • User Weight: 30 - 50 kg

#### **ADULT DIAPER (MEDIUM) (8 PCS)**

- Textile backsheet
- With barriers
- Super Absorbent paritcles
- Extra Dryness Layers
- Suitable for sensitive skin
- Special type and frontal tape system
- Wetness Indicator
- Unit per box: 8 units
- Waist Sizes: 85x125 cm
- User Weight: 40 -70 kg

[ ISO 13485 | ISO 9001

[ ISO 13485 | ISO 9001

[ ISO 13485 | ISO 9001

## **IN – PATIENT CARE**





## **Wicromed**







### Wicromed.



#### **ADULT DIAPER (MEDIUM) (30 PCS)**

- Textile backsheet
- With barriers
- Super Absorbent paritcles
- Extra Dryness Lavers
- Suitable for sensitive skin
- Special type and frontal tape system
- Wetness Indicator
- Unit per box: 30 units
- Waist Sizes: 85x125 cm • User Weight: 40 -70 kg

### **ADULT DIAPER (SMALL) (9 PCS)**

- Textile backsheet
- · With barriers
- Super Absorbent paritcles
- Extra Dryness Lavers
- · Suitable for sensitive skin
- Special type and frontal tape system
- Wetness Indicator
- Unit per box: 9 units
- Waist Sizes: 50x85 cm
- User Weight: 30 50 kg

#### **ADULT DIAPER(LARGE) (30 PCS)**

- Textile backsheet
- With barriers
- Super Absorbent paritcles
- Extra Drvness Lavers
- Suitable for sensitive skin
- Special type and frontal tape system
- Wetness Indicator
- Unit per box: 30 units
- Waist Sizes: 100x150 cm
- User Weight: 70 kg and more

#### PROTECTIVE BED SHEET

- Wicromed protective undersheet of 60 x 90 cm. rapidly absorbs wetness and transfers the same to lawyer ,hence; avoids bed to become wet and contaminated
- Whereas wicromed protective undersheet of 90 x 180 cm is designed to avoid sleep of protective undersheet during movement of the patient, and to further protect to bed
- Liquid proof PE film is pushed into edges of the bed and avoid slippage

- Unit per box: 10 unit
- Box per package: 10 boxes
- Product dimension: 60x90 cm. or 90x180 cm

- Unit per box: 30 units
- Box per package: 6 boxes
- Product dimension: 60x90 cm or 90x180cm



ISO 13485 ISO 9001

ISO 13485 | ISO 9001









#### **CLEANSING SPONGE**

- Cleansing sponge cleanses body by being bubbled by only a half glass of water
- Non allergenic
- PH value in compliance with skin (5.5 6.5)
- Coconut oil content assists regeneration of cells and removes dead cells
- Enhanced vitamins optimize skin-moisture balance
- Avoids bacteria- and fungus formation
- No rinsing is necessary
- Ideal for situations where taking a bath is difficult
- Unit per box: 20 units

#### AIR MATRESS ANTI DECUBITUS SYSTEM

#### **Pump features**

- Power supply: 220V.-240V., 50/60 H.
- **Air output:** 3-10 L.(20 -140 mmHg)
- Cycle time: 5 min
- Power consumption: 7W

#### **Mattress Features**

- Size: 186x90 cm
- 130 bubble cells
- Vents and end flaps available
- PVC thickness: 0.30 mm
- Repair kit and hose included
- Weight Capacity: 130kg
- Porduct weight: 1,5 kg

# **BLOOD PRESSURE MONITORING SYSTEM**





## **BLOOD PRESSURE MONITORING SYSTEM**



wicron (carpi



-wicron-



-wicron-PRO all.



#### **BLOOD PRESSURE MONITOR – WRIST TYPE**

- Measurement method: Oscillometric method. automatic inflation and measurement
- Blood pressure class identification
- Display: 40 x 42 mm LCD
- Date & Time display
- Measurement range Pressure: 0 300 mmHg
- Irregular cardiac rhythm dislplay (if any)
- Erroneous measurements and low battery signal
- Cuff Sizes: 300mm x 70 mm Pulse: 40 - 199 beats/min
- Accuracy Pressure: within 3 mmHg
- Pulse: within 5%
- Memories: 60 measurements memories
- **Dimensions:** 87 x 70 x 29 mm
- Weight:133 g
- Power supply: 2 x AAA alkaline batteries
- Exclusice storage case
- Operating temperature +10 to +40°C, 85% R.H.
- Storage temperature -20 to +50°C, 85% R.H.



**BLOOD PRESSURE MONITOR – WRIST TYPE** 

- Measurement method: Oscillometric method. automatic inflation and measurement
- Blood pressure class identification Guiding Human voice . 3 language options
- Display: 40 x 42 mm LCD
- Date & Time display
- Measurement range Pressure: 0 300 mmHg
- Irregular cardiac rhythm dislplay (if any)
- Erroneous measurements and low battery signal
- Cuff Sizes: 300mm x 70 mm • Pulse: 40 - 199 beats/min
- Accuracy Pressure: within 3 mmHg
- Pulse: within 5%
- **Memories:** 60 measurements memories
- **Dimensions:** 87 x 70 x 29 mm
- Weight: 133 g
- **Power supply:** 2 x AAA alkaline batteries
- Exclusice storage case
- Operating temperature +10 to +40°C, 85% R.H.
- Storage temperature -20 to +50°C, 85% R.H.
- Talking feature( Arabic Kurdish Russian )

ISO 13485 | ISO 9001

#### **BLOOD PRESSURE MONITOR – ARM TYPE**

- Measurement method: Oscillometric method. automatic inflation and measurement
- Blodd pressure class identification
   Guiding Human voice , 3 language options
   Display: 70 x 75 mm LCD
- Date & Time display
- Measurement range Pressure: 0 295 mmHg Pulse: 40 199 beats/min
- Irregular cardiac rhythm dislplay (if any)
- Cuff circumference:

22cm-30cm(8 21/32" -11 13/16")
30cm-42cm(11 13/16" -16 17/32") (Optional)
42cm-48cm(16 17/32" -18 29/32") (Optional)
• Erroneous measurements and low battery signal
• Cuff Sizes: 300mm x 70 mm

- Accuracy Pressure: within 3 mmHg
- Pulse: within + 5%
- Memories: 60 measurements memories
- Dimensions 125 x 130 x 61 mm **Weight:** 330 gr
- Power supply: 4 x AA alkaline batteries
- Exclusice storage case
- Operating temperature +10 to +40°C, 85% R.H. Storage temperature -20 to +50°C, 85% R.H.

### -wicron-Pump



### -wicron-PRO



### -wicron-nurse



#### **BLOOD PRESSURE MONITOR – ARM TYPE**

- 48 mm. scale
- 0-300 mm. /Hg. dial
- Membrane protected single hose output
- 3 mm. /Hg. accuracy
- Stethoscope
- Adhesive or hook cuff option
- Leather storage bag
- 2 year warranty

#### **BLOOD PRESSURE MONITOR FULLY** AUTOMATIC - ARM TYPE

- Measurement method : Oscillometric method, automatic inflation and measurement
- Blodd pressure class identification
   Display 70 x 75 mm LCD

- Date & Time display
  Measurement range Pressure: 0 295 mmHg Pulse: 40 - 199 beats/min
- Irregular cardiac rhythm dislplay (if any)
- Cuff circumference:

22cm-30cm(8 21/32" -11 13/16")

30cm-42cm(11 13/16" -16 17/32") (Optional) 42cm-48cm(16 17/32" -18 29/32") (Optional)

- Erroneous measurements and low battery signal
- Cuff Sizes: 300mm x 70 mm
- Accuracy Pressure: within 3 mmHg
   Pulse: within ± 5%
- Memories : 60 measurements memories
- Dimensions 125 x 130 x 61 mm
- Weight :330 gr
- Power supply: 4 x AA alkaline batteries
- Exclusice storage case
- Operating temperature +10 to +40°C, 85% R.H.
- Storage temperature -20 to +50°C, 85% R.H.

ISO 13485 | ISO 9001

#### MANUAL BLOOD PRESSURE MONITOR WITH **CUFF - NURSE TYPE**

- · Gear-free trigger aneroid
- Certified accuracy to ±3 mmHg
- Durable: withstands a 30" drop
- · Laser-engraved dial face for accuracy
- Measrument range : 0 − 300 mmHg
- Black PVC hose
- Soft comfort grip
- Trigger deflation valve for smooth, one-finger operation
- Cuff tested to withstand demanding clinical environments
- Rotatable port, improves patient comfort and cuff durability
- · Folded cuff edge reduces risk of cuts and scrapes
- Cuff meets all clinical guidelines for proper fit from AAMI and AHA
- Dimension: 73mm x 71 mm x 84 mm
- Weight: 272 gr
- Exclusive case (Black)

ISO 13485 | ISO 9001







## Wicromed.



**CERVICAL COLLAR - FOAM** 

#### Indications

- Acute neck pain
- Cervical disc hernia
- Minor muscle spasm associated with spondylosis
- Protection during halo application

#### **Functions**

- Kinesthetic reminder to limit neck movement
- · Limitation of movement

#### **Features**

- · Made of foam with a stockinette cover
- · High density foam is available



## Wicromed.



**CERVICAL COLLAR - NELSON** 

#### Indications

- Cervical disc hernia
- Post-traumatic protective purpose and patient transfer
- Cervical strain
- Postural stiffness

#### **Functions**

· Vents which ensure air circulation

ISO 13485 | ISO 9001

Kinesthetic reminder to limit neck movement

Made of antiallergenic 12 mm. Plastoze material

### Wicromed.



**CERVICAL COLLAR - PHILADELPHIA** 

#### Indications

- Post-operative cervical fractures
- Halo removal
- Anterior cervical fusion
- Dens type I cervical fracture of C2
- Anterior discectomy
- Suspected cervical trauma in unconscious patients
- Tear-drop fracture of the vertebral body (Note: Some tear-drop fractures require anterior decompression and fusion.)
- Cervical strain

#### Functions

- · Kinesthetic reminder to limit neck movement
- Plastic struts restrict neck flexion, extension, lateral bending and rotation

- Made of antiallergenic 12 mm. Plastoze material
- An anterior hole for a tracheostomy is available





#### **CLAVICULA BANDAGE**

#### Indications

- Poor posture
- Clavicula fracture
- Clavicola dislocation
- Sterno-clavicular dislocations
- · Weak scapula muscles

#### Functions

- Excellent aid in correcting poor posture and shoulder slump
- Help to correct and remind poor shoulder posture
- Immobilize fractured clavicle

#### **Features**

- Made of antiallergenic sponge coated by cotton fabrics
- Standard

### ISO 13485 | ISO 9001

## Wicromed.



#### **SHOULDER SUPPORT (NEOPREN)**

#### Indications

- Shoulder pain
- Bursitis
- Tendinitis
- Shoulder articulation semi dislocation

## ISO 13485 | ISO 9001

## Wicromed



#### **VELPEAU SHOULDER SUPPORT**

#### Indications

- S/P glenohumeral acromioplasty
- Surgical repair of the rotator cuff
- Anterior glenohumeral capsular reconstruction
- Anterior dislocation of the G-H joint
- Swelling and inflammation
- Torn muscles or ligaments
- Impingement syndrome
- Tendonitis
- Myositis

#### **Functions**

- Protects your arm from rotating
- Restricts mobility of the shoulder and upper arm for separations, dislocations, strains and other injuries

#### Features

- Breathable cotton material
- Adjustable hook and loop closures allow swathe to be positioned in desired location
- Universal design fits left and right arm



ISO 13485 | ISO 9001

## Wicromed.



### **ARM SLING (30 DEGREE)**

#### Indications

- The impingement syndrome of post-op
- Rotator cuff tear after repair.
- Frozen shoulder under anesthesia after manipulation.
- Subacromial bursitis and subkorakoid.
- Used after shoulder endoprosthesis.
- Bursitis

#### Functions

Abduction of shoulder articulation at 30 degrees

#### Features

- Inner surface cotton towel, outer surface Velcro: made of sponge laminated by maintaining fabrics
- Adjustable tab



## Wicromed



#### **ARM SLING - GABARDINE**

#### Indications

- Brachial nerve injury
- Upper arm fracture
- Shoulder subluxation/ dislocation
- Post surgical protection

#### **Features**

- Made of cotton fabricsSupports the arm and carries the weight across the back and shoulder
- Comfortable lightweight material
- · Immobilization of shoulder

## Wicromed.



#### **ARM SLING -WITH NET**

#### Indications

 Dislocations in shoulder and half of the sholder dislocations. In cases of postoperative resting, Fractures of the distal edema formation is used to prevent.

#### Features

- . The inner surface of the towel, the exterior is made of airtex cloth. The waist belt provides extra protection supportes.
- Black

## Wicromed.



#### ARM SLING -SPONGE

#### Indications

- Brachial nerve injury
- Upper arm fracture
- Shoulder subluxation/ dislocation
- Post surgical protection

#### Features

- Inner surface cotton towel, outer surface Velcro: made of sponge laminated by maintaining fabrics
- Supports the arm and carries the weight across the back and shoulder
- Universal design fits right or left arm.
- Adjustable slide buckle for proper height adjustment
- · Comfortable lightweight material

### Wicromed.



#### **ELBOW BRACE - EPICONDYLITIS** SUPPORTED

#### Features

 Circular pad for extra pressure on the muscle in the arm or wrist

#### Indications

- Tennis elbow
- · Golfer's elbow
- · Weak elbow
- · Articular effusion & swelling
- Chronic inflammation after operation or trauma
- · Collateral ligament injury

#### **Functions**

- Closed cell neoprene with compressive and proprioceptive effect
- Increase blood circulation through warm effect









#### **Features and Functions**

**EPICONDYLITIS BANDAGE** 

- Silicone pad provides compression to help relieve tennis elbow pain and discomfort
- Helps the body quickly reabsorb edema and hematoma
- Relieves strain on the periarticular tissues and compresses the soft tissue

#### Indications

- · Tennis elbow
- · Golfer's elbow
- Elbow traumas
- · Ligament injuries, post-operation during healing, Sport activities.

## Wicromed.



#### **EPICONDYLITIS BANDAGE (DELUX)**

- · Used in epicondylitis treatment
- Decrease pressure on tendon by compression

## Wicromed.



#### **DEQUARVEIN SPLINT**

#### **Features and Functions**

- Removable metal splints in thumb area for improved
- Supports the wrist and thumb, leaving the palm free
- Ideal for sprains of thumb, ligament support, and minor fractures of the thumb

## Wicromed.



#### HAND WRIST SPLINT

#### Indications

- Carpal Tunnel Syndrome
- Wrist instability due to triangular fibrocartilage complex TFCC injury
- Radial nerve paralysis
- Rheumatoid arthritis

#### **Functions**

- Immobilization of the wrist
- · Firm support for wrist

#### Features

- · Comfortable, breathable, skin-friendly flannel-lined elastic cotton fabrics
- Anatomically-contoured aluminium palmar splint

[ ISO 13485 | ISO 9001

ISO 13485 | ISO 9001

[ ISO 13485 | ISO 9001

## Wicromed.



#### HAND WRIST SPLINT - AIRTEX

#### Indications

- In the case of carpal tunnel and along dequervain
- · Dequervain tenosybovitis
- Red Blue Black

#### Feeatures and Functions

- Enables immobilization of the wrist in pain
- Inhibits movement of the wrist with the palmar aluminium
- · Made of airtex fabric.
- · Extremely breathable
- · Lightweight, durable soft fabric
- Anatomically contoured aluminum supporter for wrist
- · Available in right and left hand

## Wicromed.



HAND WRIST SPLINT WITH THUMB SUPPORT AIRTEX

#### Indications

- Carpal Tunnel Syndrome
- De Quervain's disease
- First metacarpophalangeal or carpometacarpal joint injury
- DeQuervain's tenosynovitis & Gamekeeper's thumb
- Ligamentous laxity of the thumb skier's thumb
- Radial nerve paralysis
- Rheumatoid arthritis or degenerative arthritis at the thumb
- Red Blue Black

#### **Functions**

Immobilization of the wrist and thumb

#### Features

- Breathable material with terry-lined interior for wearing comfort
- Made of elastic cotton fabric
- Anatomically-contoured aluminium palmar splint
- Rigid thumb stay for custom fit
- V-notch at fifth MP for palmar spread
- · Available in right and left hand



## Wicromed.



#### HAND WRIST SPLINT - NEOPREN

#### Indications

- Ideal for sprains and strains and other wrist injuries.
   Features and Functions
- Wrist immobilization
- · Decrease anti-inflammatory medications
- Decrease pain symptoms
- Improve wrist joint function
- Easy application

ISO 13485 | ISO 9001

• Inhibits movement of the wrist by Palmar aluminum support.

### Wicromed.



HAND WRIST SPLINT WITH THUMB SUPPORT

#### Indications

- Carpal Tunnel Syndrome
- First metacarpophalangeal or carpometacarpal joint injury
- DeQuervain's tenosynovitis & Gamekeeper's thumb
- Ligamentous laxity of the thumb skier's thumb
- Radial nerve paralysis
- Rheumatoid arthritis or degenerative arthritis at the thumb

#### Functions

· Immobilization of the wrist and thumb

#### Features

- Breathable material with terry-lined interior for wearing comfort
- Anatomically-contoured aluminium palmar splint
- · Rigid thumb stay for custom fit
- V-notch at fifth MP for palmar spread
- · Available in right and left hand





#### HAND WRIST SPLINT WITH THUMB SUPPORT - NEOPREN

#### Indications

- Carpal Tunnel Syndrome
- First metacarpophalangeal or carpometacarpal joint injury
- DeQuervain's tenosynovitis & Gamekeeper's thumb
- Ligamentous laxity of the thumb skier's thumb
- Radial nerve paralysis
- Rheumatoid arthritis or degenerative arthritis at the thumb

#### **Functions**

· Immobilization of the wrist and thumb

#### Features

- Breathable material with terry-lined interior for wearing comfort
- Anatomically-contoured aluminium palmar splint
- · Rigid thumb stay for custom fit
- V-notch at fifth MP for palmar spread
- · Available in right and left hand

**C** € ISO 13485 | ISO 9001

## Wicromed.



#### **WOOLEN CORSET**

#### Indications

· Rheumatism, sciatica, lumbago

#### **Features and Functions**

Made of elastic wool fabrics

ISO 13485 ISO 9001

- Made of wool fabrics 30 cm high
- To maintain warm both the waist and abdominal muscles

# Wicromed.



ABDOMINAL SUPPORT CORSET ( HOOKED)

#### Indications

- After birth
- Mild low back pain
- Splanchnoptosis and exomphalos

#### **Features and Functions**

- Provides strong comfortable support to help to relieve back pain during pregnancy
- Specially-shaped to provide full support in the sacral region and over the symphysispubis

## Wicromed.



#### **CHEST CORSET (HOOKED)**

#### Indications

- To support sternum post-coronel surgery
- Post-Costa fracture

#### **Features and Functions**

- Made of air permeable, antiperspirant fabrics with nylon thread
- Additionally supported by towel laminated soft sponge fabrics on sternal site

ISO 13485 | ISO 9001

## Wicromed.



#### **CORSET LUMBOSTAD (HOOKED)**

#### Indications

- Lumbago
- Spinal deformations
- Sciatica pain
- Lumbal vertebral problems
- Acute or chronic lumbago
- Post-up immobilization purposed
- · Thermo-therapy purposed

#### **Features and Functions**

- Made of antiallergenic, ventilated, antiperspirant, elastic cotton material: made of 26 cm in width
- Lumbal site is supported by plastic and steel underwires

### Wicromed.



**CORSET LUMBO-SACRAL (HOOKED)** 

#### Indications

- Lumbago
- Spinal deformations
- Vertebral and ilio-sacral articulation problems
- Acute or chronic lumbago
- Post-up immobilization purposed

#### **Features and Functions**

- Made of antiallergenic, ventilated, antiperspirant, elastic cotton material; made of 32 cm in width,
- Lumbal and sacral sites are supported by plastic and steel underwires,
- Supported by additional belts for more compliance with waist
- Plastic underwire 2 at front, 4 at back
- . Optionally; 2 steel underwires may be added to underwires at back, which ensures fully stabile structure

### ISO 13485 | ISO 9001

## Wicromed.



#### **CORSET LUMBO-SACRAL**

#### Indications

- Lumbago
- Spinal deformations
- Vertebral and ilio-sacral articulation problems
- · Acute or chronic lumbago
- Post-up immobilization purposed

#### **Features and Functions**

- Made of antiallergenic, ventilated, antiperspirant, elastic cotton material; made of 32 cm in width,
- Lumbal and sacral sites are supported by plastic and
- Supported by additional belts for more compliance with waist.
- Plastic underwire 2 at front, 4 at back
- . Optionally; 2 steel underwires may be added to underwires at back, which ensures fully stabile structure

## ISO 13485 | ISO 9001

## Wicromed.



#### **CORSET POSTUREX (HOOKED)**

#### Indications

- Poor posture
- Weak scapula muscles
- Round shoulder

#### Features and Functions

- Triangle configuration correctly distributes pressure to tip of shoulder
- Additional stays help reinforce back alignment and proper posture
- · Crisscrossing straps provide more support and make the support easy to adjust
- Corrects poor posture and prevents slouching







#### **CORSET DORSOLOMBER (HOOKED)**

#### Indications

- Lumbargo
- Osteoporosis
- · Light weight
- Post-up immobilization purposed

#### Features and Functions

- Made of antiallergenic, ventilated, antiperspirant fabrics with nylon thread
- Supported by thermoplastic and steel underwires on thoroca, lumbal and sacral sites
- Supported by additional belts for more compliance with waist and breast
- Plastic underwires 2 at front, 4 at back
- Optionally; 2 steel underwires may be added to underwires at back, which ensures fully stabile structure

## Wicromed.



#### **HERNIA TRUSS – SINGLE SIDED**

#### Indications

- Mild lumbago
- Sciatica pain
- · Bilateral inquinal hernia
- Mild scoliosis

#### Features and Functions

- Provides compression around the affected area with targeted compression on the hernia, reducing the inguinal hernia
- Hook and loop closure allows to controlled fits
- Adjustable thigh strap are made of soft elastic

## Wicromed



#### HERNIA TRUSS – DOUBLE SIDED

#### Indications

- Bilateral inquinal hernia
- Mild lumbago
- Sciatica pain
- Mild scoliosis

#### Features and Functions

- Provides compression around the affected area with targeted compression on the hernia, reducing the inguinal hernia
- Infinite adjustment for maximum comfort minimizes pressure over back musculature
- A simple hook and loop closure allows you to control the fit and compression of the truss
- Elastic band trusses should be put around the waist with the pad places in position and tightened to fit

### Wicromed.



### **CALF SUPPORT (UP)**

#### Indications

- Muscle weakness
- Muscle fatigue

#### **Features and Functions**

- Helps relieve pain from an injured or weak thigh
- Provides heat retention
- Effective support for strained groin muscles and torn hamstring



ISO 13485 ISO 9001

 $C \in$ 

ISO 13485 | ISO 9001



## Wicromed.



#### **CALF SUPPORT (UP)**

#### Indications

- Muscle weakness
- Muscle fatique or strain

#### **Features and Functions**

- Helps relieve pain from an injured or weak thigh
- Provides heat retention
- Effective support for strained groin muscles and torn hamstring

## Wicromed.



#### **CALF SUPPORT (DOWN)**

#### Indications

- Muscle weakness
- Compartment syndrome
- Thermal therapy purposes

#### Features and Functions

- Helps relieve pain from an injured or weak shin
- Designed to conform naturally to the contours of the shin
- Effective support for strained muscles and torn hamstrings

## Wicromed.



#### **CALF SUPPORT (DOWN)**

#### Indications

- Tibial bones protection case
- Muscle weakness
- Compartment syndrome
- Thermal therapy purposes

#### **Features and Functions**

- Provide warm and compression in muscular lesion
- Helps relieve pain from an injured or weak shin

### Wicromed.



#### **WOOLEN KNEE SUPPORT**

#### Indications

- Thermal-therapy purpose
- Knee compression
- Arthritis
- Mild strain/spain

#### **Features and Functions**

- Made of wool fabrics of 30 cm high
- For knee pains, rheumatism and sciatic, and where knees must be maintained warm
- Unnoticeable under clothes
- Does not irritate skin and help prevent allergies







ISO 13485 | ISO 9001





#### **KNEE SUPPORT - NEOPRENE**

#### Indications

- Knee trauma
- Osteoarthritis
- Articular effusion
- Local inflammation
- Thermo-therapy purpose

#### **Features and Functions**

- Helps relieve weak or sore knees
- · Side air vents for fit and comfort.
- Comfortable neoprene blend helps knee retain warmth.
- Adjustable straps for proper support, comfort and fit.
- Easy to apply

### [ ISO 13485 | ISO 9001

## Wicromed.



#### **CLOSED PATELLA KNEE SUPPORT – NFOPREN**

#### Indications

- Articular effusion & swelling
- Chronic inflammation after operation or trauma
- Patellofemoral pain syndrome
- Bilateral spiral steel stays prevents valgus or varus deviation
- Patella alta, patella subluxation, patella dislocation
- Chondromalacia
- · Muscle strain around knee joint

#### **Functions**

- · Effectively stabilizes the kneecap to eliminate anterior
- Closed cell neoprene with compressive and propriocep-
- Increase blood circulation through warm effect

#### **Features**

- 1/8" thick neoprene with nylon lining
- Open patellar design
- Popliteal cut-out design



ISO 13485 | ISO 9001

## Wicromed



#### PATELLA TENDON STRAP

#### Indications

- Patellar tendinopathy
- Osgood-Schlatter disease
- Jumper's knee

#### **Features and Functions**

[ ISO 13485 | ISO 9001

- Reinforces and protects the patellar tendon
- 1/8" thick neoprene with nylon lining
- Simple fitting & time-saving application via hook and loop closure

### Wicromed.



#### PATELLA KNEE SUPPORT

#### Indications

- Osteoarthritis
- Articular effusion & Swelling
- Femoral pain

#### Features and Functions

- · Stabilize patella by special ped support
- Effectively stabilizes the kneecap to eliminate anterior
- Increase blood circulation through warm effect
- Simple fitting & time-saving application via hook and loop closure



## Wicromed.



#### PATELLA KNEE SUPPORT

#### Indications

- Articular effusion & swelling
- Chronic inflammation after operation or trauma
- Patellofemoral pain syndrome
- Patella alta, patella subluxation, patella dislocation
- Chondromalacia
- · Muscle strain around knee joint

#### **Functions**

- Effectively stabilizes the kneecap to eliminate anterior
- Closed cell neoprene with compressive and proprioceptive effect
- · Increase blood circulation through warm effect

#### Features

- 1/8" thick neoprene with nylon lining
- Open patellar design
- · Popliteal cut-out design



## $\mathcal{W}$ icromed



#### PATELLA AND LIGAMENT KNEE SUPPORT (LONG - SHORT)

#### Indications

- Mild instabilities of knee joint
- Articular effusion & swelling
- Chronic inflammation after operation or trauma
- Patellofemoral pain syndrome
- Patella alta, patella subluxation, patella dislocation
- Chondromalacia
- · Muscle strain around knee joint

#### Functions

- Bilateral spiral steel stays prevents valgus or varus
- Effectively stabilizes the kneecap to eliminate anterior knee pain
- Closed cell neoprene with compressive and proprioceptive effect
- Increase blood circulation through warm effect

#### **Features**

- 1/8" thick neoprene with nylon lining
  Crossing straps control patella movement.
- Open patellar design with soft-cushioned felt
- Upper and lower coronal fasteners offer further fit
- Pull-on version

ISO 13485 | ISO 9001

## Wicromed



#### PATELLA AND LIGAMENT KNEE SUPPORT

#### Indications:

- Arthritis
- Mild instabilities of knee joint
- · Articular effusion & swelling
- Chronic inflammation after operation or trauma
- Patellofemoral pain syndrome
- Patella alta, patella subluxation, patella dislocation
- Chondromalacia
- · Muscle strain around knee joint

#### Features and Functions

- Neoprene circle pad around the patella provides padded protection and stabilization to prevent patella displacement
- Two medial and lateral stays provide additional support and stability for the knee
- Helps relieve pain from tendinitis and arthritis
- Increase blood circulation through warm effect



## Wicromed



#### HINGED PATELLA KNEE SUPPORT

#### Indications

- · Mild to moderate instabilities of knee joint
- Mild to moderate strain of the collateral ligaments
- Meniscus iniuries
- Chronic inflammation after operation or trauma

- Collateral stabilization of knee joint prevents valgus or varus deviation.
- Prevent hyperextension
- Closed cell neoprene with compressive and propriocentive effect
- Increase blood circulation through warm effect

#### **Features**

- Open patellar design
- Hinges hidden in pockets offer superior stability
- Soft-cushioned felt forms the contour of knee and offers additional comfort.
- Upper and lower coronal fasteners offer further fit.
- Pull-on version





#### **HINGED PATELLA KNEE SUPPORT (OPEN** FRONT)

#### Indications

- Weak knee
- Moderate arthritis
- Chronic inflammation after operation or trauma
- Patellofemoral pain syndrome
- Patella alta, patella subluxation, patella dislocation
- Knee swelling

#### Features and Functions

- The front opening design and easy open closure straps make the brace suitable of the arthritic or elderly patient who may suffer from limited dexterity
- Elastic brace with motion hinges provides medial and lateral support for the unstable knee
- May be used for stabilization and support of the knee following strains, injury or surgery
- Open patellar design
- Popliteal cut-out design



## Wicromed.



#### HINGED KNEE BRACE

#### Indications

- Mild to moderate instabilities of knee joint
- Mild to moderate strain of the collateral ligaments
- Meniscus iniuries
- Chronic inflammation after operation or trauma

#### **Functions**

- Collateral stabilization of knee joint prevents valgus or varus deviation
- Prevent hyperextension
- Closed cell neoprene with compressive and proprioceptive effect
- Increase blood circulation through warm effect

- 3/16" thick neoprene with nylon lining
- Open patellar design
- Two hook and loop straps insure comfortable fit
   Soft-cushioned felt forms the contour of knee and offers additional comfort
- Upper and lower coronal fasteners offer further fit.
- Pull-on version
- Provides firm support for the maintenance of joint structures with gentle compression over the kneecap



## Wicromed.



#### ANKLE SUPPORT WITH STRAP

#### Indications

- Mild to moderate ankle distortions
- Treatment of acute ankle sprain (alternative to taping)
- · Chronic ankle instability
- Articular effusion & swelling
- Prophylactic and therapeutic use in sport activities

#### **Functions**

- Prevent over inversion or plantarflexion
- Closed cell neoprene increases blood circulation through warm effect
- · Compressive and proprioceptive effect

#### Features

 Extra elastic figure 8 strap provides additional compression and support

### ISO 13485 | ISO 9001

## Wicromed.



#### **ANKLE SUPPORT**

#### Indications

- · Mild torsion of ankle joint
- Articular swelling and chronic ankle instability

#### **Functions**

- Fits well under clothing
- Helps relieve pain due to strains or stiffness
- Compressive and proprioceptive effect

#### **Features**

· Can also for prophylactic use in low intensity exercise or general workout



## Wicromed.



#### ANKLE SUPPORT - EXTRA

#### Indications

- Thendomyopathies
- Arthritis-arthritis
- · Situations where stabilization is necessary
- Thermal-therapy purposed.

#### **Functions**

- Designed to protect ankle Malleolus
- 2 asymmetric designs formed in anatomical structure
- Compresses the site without avoiding motion
- Structure which ensures rapid absorption of blood circulation and edema, discharge and hematoma

## Wicromed.



#### ANKLE SUPPORT WITH SILICONE PAD

#### Indications

- Mild ankle distortion
- Chronic ankle instability
- Articular effusion & swelling
- Prophylactic and therapeutic use in sport activities

#### **Functions**

- The silicone insert exert a massage effect.
- Reduces swelling from minor injuries
- Flexible compression provides ankle support and improves proprioception

#### **Features**

- Silicone support pad
- Made of seamless, four-way stretch elastic material for comfortable fit
- Breathable material with moisture absorbed



## Wicromed.



#### ANKLE SUPPORT WITH LIGAMENT

- Provides compression for weak or overstressed ankles
- Thermal-therapy purposed.

### Wicromed.



#### AIRCAST ANKLE SUPPORT

#### Indications

- · Mild to moderate ankle sprain
- Chronic ankle instability
- Prophylactic and therapeutic use in sport activities

#### **Functions**

- Rigid contoured outer shells control inversion or eversion while allowing for normal dorsi-flexion
- Inflatable air bladders allow user to adjust compression pressure to reduce swelling

#### Features

- Air bladders for additional support and compression
- Anatomically-contoured plastic shell
- Simple to apply and adjust by two circumferential
- Easy adjustable heel strap
- · Designed specifically for right or left



ISO 13485 | ISO 9001





## Wicromed.



#### **HALLUX VALGUS BANDAGE**

#### Indications

- Hallux valgus
- Features and functions
- Aluminum underwire on lateral site of thumb
- Standard (Right-Left)

#### STRESS BALL

- Made of 100% silicone, different color and hardness
- Rehabilitation oval silicone ball is used for exercise.
- It has got 3 grade red (hard), blue (medium), green (soft)

# **METAL ORTHOPEDICS**





# **METAL ORTHOPEDICS**



### Wicromed.



# Wicromed.



### Wicromed.



### **WHEELCHAIR MANUEL - WM 809**

• Folding Width: 23 cm Seat width: 46 cm • Rear Wheel Dia: 59 cm • Lower Knee Length: 15 - 51 cm • Seat depth: 40 cm • Back height: 40 cm • Total Load: 100 KG • Weight: 19 Kg

• Other: leather designed , blue

• Chair Width: 60 cm

**( ( ISO** 13485 **| ISO** 9001

### POWER WHEELCHAIR - HBLD2 -A

Powder coated steel folding frame, Flip up armrest, detachable footrest, drop back handle, 200x50 PU castors

- Reinforced aluminium body ,readily portable by folding option; double motor wheelchair
- Electronic Joystick ensures comfortable and safe drive
- Weight (battery excluded): 24 kg • Dimensions (LxWxH): 110x64x91 cm
- Weight (battery include ): 51 kg
- Maximum load capacity: 100-120 kg
- Max. Speed: 6 km/h
- Escalation angle: 10 degrees
- Serially connected 2 batteries: 2x12 V. 20 Ah
- Motor power: 2 x 220 W
- Operation time (full charge): 4 hours
- Noise level: 60 db

### WHEEL CHAIR - WM 607-GCJ

Chromed steel frame, 8" PVC castor, 24" Solid wheel. detachable armrest, double control brake system, adjustable footrest, Reclining high back, chamber pot.

- Chair Width: 65 cm
- Folding Width: 30 cm • Seat width: 46 cm
- Rear Wheel Dia: 59 cm
- Lower Knee Length: 55 cm
- Seat depth: 43 cm • Back height: 74 cm
- Total Load: 90 Kg • Weight: 31 Kg

[ ISO 13485 | ISO 9001



ISO 13485 ISO 9001

# Wicromed.



### WHEELCHAIR - TRANSPORT- HOME TYPE WM 976 -AJ -43

• Seat width: 43 cm • Seat depth: 40 cm • Back height: 54 cm • Total Load: 100 Kg • Folding Width: 28 cm • Weight: 22 Kg

• Lower Knee Length: 51 cm • Rear Wheel Dia: 20 cm

 Other: Chrome coated , steel frame , fixed armrest. removable footrest, double control braking system

# Wicromed.



### WHEELCHAIR - FOOTREST RAISABLE WM 902 -C

• Chair Width: 53 cm • Folding Width: 28 cm • Seat width: 41 cm • Rear Wheel Dia: 59 cm • Lower Knee Lenath: 15-50 cm

• Seat depth: 40 cm

• Back height: 54 cm • Total Load: 100 Kg • Weight: 22 Kg

• Other: Black leather, metalic silver color, footrest raisable, armrest removable

# Wicromed.



### **COMMODE WHEELCHAIR - WITH TOILET CHAMBER WM 608**

- Chair Width: 65 cm • Folding Width: 24 cm
- Seat width: 44 cm
- Rear Wheel Dia: 59 cm
- Lower Knee Lenath: 52 cm
- Seat depth: 45 cm • Back height: 45 cm
- Total Load: 100 Kg
- Weight: 20 Kg
- Other: Chrome coated, steel frame, fixed armrest and footrest

### Wicromed.



### MANUEL WHEELCHAIR FOR CHILDREN WM 981L - 35

- Chair Width: 35 cm
- Folding Width: 25 cm • Seat width: 38 -55 cm
- Rear Wheel Dia: 56 cm
- Lower Knee Length: 45 cm
- Seat depth: 36 cm
- Back height:38-55 cm
- Total Load: 75 Kg
- Weight: 14 Kg
- Other: Aluminum coating; for children use; light components, removable footrest, adjustable backrest.

[ ISO 13485 | ISO 9001

[ ISO 13485 | ISO 9001

[ ISO 13485 | ISO 9001

**( ( I I SO 13485 | I SO 9001** 

# **METAL ORTHOPEDICS**

# Wicromed.



### MANUEL WHEELCHAIR FOR CHILDREN WM 802 - 35

- Width: 48 cm • **Depth:** 58 cm • **Height:** 76 – 85 cm • Weight: 5 Kg
- Other: Knob . metalic silver color . fixed-portable . foldable

# Wicromed.



### MANUEL COMMODE WHEELCHAIR - WITH **NECK STABILIZER 609 -GC**

- Chair Width: 65 cm • Folding Width: 33 cm
- Seat width: 46 cm • Rear Wheel Dia: 59 cm
- Lower Knee Lenath: 15 54 cm
- Seat depth: 43 cm • Back height:74 cm • Total Load: 90 Kg
- Weight: 14 Kg
- Other: Black leather designed, metalic silver color, detechable commode . Neck stabilizer

### Wicromed.



### MANUEL WHEELCHAIR - WIDER TYPE -WM 950 LBQ 46

- Chair Width: 65 cm
- Folding Width: 34 cm
- Seat width: 46 cm
- Rear Wheel Dia: 61 cm
- Lower Knee Length: 52 cm
- Seat depth: 40 cm
- Back height: 40 cm
- Total Load: 100 Kg
- Weight: 16 Kg
- Other: Aluminium coated , soft cushion ,8" PVC castor removable footrest, adjustabl brake system

### Wicromed.



### CRUTCH - WITH SCREW

- Soft rubber and knob
- Adjustable height
- Modern design

	Weight	Height
Small	1,10 kg	90 - 105 cm
Medium	1,30 kg	110 -130 cm
Larne	1.50 kg	120 -140 cm

[ ISO 13485 | ISO 9001

[ ISO 13485 | ISO 9001

[ [50 13485 | 150 9001

( F ISO 13485 | ISO 9001

# Wicromed. **CRUTCH - WITH CLIP** Soft rubber and knob Adjustable height



• Weight: 0,80 kg

Wicromed.

• Heght: 97 -118

FOREARM CRUTCH

• Other: Adjustable height, 6 color options, knob,

# **Wicromed**.



Other: Metalic silver color , movable , foldable , shiny

# Wicromed.



#### **WALKER - ALUMINIUM WALKER - METAL**

- Width: 53 cm • Width: 53 cm • **Height:** 76 -85 cm • **Height:** 76 -85 cm • **Dept:** 40 cm • **Dept:** 40 cm • Weight: 3,10 kg • Weight: 3,10 kg
  - Other: Metalic silver color, movable, foldable, shiny anodizing aluminium

[ ISO 13485 | ISO 9001

Modern design

Small

Large

Medium

Weiaht

1,40 kg

1,50 kg

1,60 kg

Height

90 - 105 cm

110 -130 cm

120 -140 cm


ISO 13485 | ISO 9001


[ ISO 13485 | ISO 9001

anodizing aluminium

[ ISO 13485 | ISO 9001

# **METAL ORTHOPEDICS**





# **METAL ORTHOPEDICS**

# Wicromed.



### **UPLIFT COMMODE ASSIST**

- Adaptable with all type of toilet
- Cover or without cover options are available
- Uplift the seat height for 11 cm
- Special clips prevent slips

# **FOOT HEALTH**





# **FOOT HEALTH**



Wicromed.



Wicromed.



Wicromed.



### SILICONE ORTHOPEDIC INSOLE

- Made of visco elastic silicone
- Hypo-Allergenic
- Relaxing and soothing
- Ensures balanced pressure distribution in foot sole by supporting the same
- anatomically
- Blue regions are softer compared to others
- Extends until metatarsus
- % 100 silicone
- Elastic structure

### SILICONE INSOLE (HIGH ARC SUPPORTED)

- · Made of 100% silicone
- Absorbs ground reaction force
- Orthopedic supports ensure balanced load distribution Indications:
- Pet planus
- Pes caus
- Plantar fasilit
- Metetarsal support

### **CORK INSOLES**

For lifting the medial longitudinal arch;uplifting of the back part of the root,adduction of the forefoot. Custom made on cast by molding high quality leather and cork. Indications:

- Straight sole
- Pes planus
- Pes cavus
- Metatarsalji

Size Range 33-34

**( €** ISO 13485 | ISO 9001

**( €** | ISO 13485 | ISO 9001

**C** € 150 13485 | 150 9001

### Wicromed. **Wicromed**. Wicromed. Wicromed. SILICONE HEEL - FABRIC COVERED STEEL INSOLE SILICONE EPIN HEEL SILICONE METATARSAL PAD (ADHERENT) • Supported by stainless steel bridge; features appearent Made of 100% silicone . Made of visco elastic silicone Made of 100% silicone Hypo-Allergenic Helps reduce and absorb the shock with each step. Indications Maintains extended term elaticity and Dynamics thereof Blue regions are softer compared to others Ideal for healing pains Metetars head sensitivity - Metatarsalgia Made of naturel leather Epin calcanei For epin calcai Fat Pad Atrophy One size smaller than the foot must be used Lower-Metetars callus % 100 silicone • Ideal for conditions such as Fat Pad Atrophy Indications: Washable and reusable - Easily fits into all shoes. Pes planus Pes cavus Metatarsalji Elasticity Size Range 34-36

[ ISO 13485 | ISO 9001

# **FOOT HEALTH**

### Wicromed




### Wicromed.



### Wicromed.



### MEDI **PHOENIX**



### SILICONE BUNION CUSHION

- Made of 100% silicone
- Helps with bunion discomfort, pressure and friction
- Helps protect/cushion protuberance
- Premium Quality Silicone for comfort and durability
- Ergonomically designed to fit the anatomy of the wearer's

#### SILICONE TOE SPREADER

- . Made of visco-elastic silicone
- Designed slimmer in the middle for greater comfort, broader top and bottom to maintain spread
- Keeps toes apart and reduces friction
- Reusable.

#### Indications

- Toe deformities
- · Reduces friction & irritation.
- · Provides gentle, constant pressure that promotes proper joint alignment.
- Aligns the big toe to help relieve bunions

### **MEDICAL CLOGS**

- · Perforated full anatomic Clogs
- Hygienic and does not produce bacteria
- Availabel for blue and White colour.

MAN	WOMAN
40 -45	35 -39

### **VARICOSE SOCKS ANTIEMBOLIC UPPER / LOWER KNEE**

- Prevention of Pre/post-operation thrombosis and
- Situations where general immobility, bedridden
- Pregnancy, birth and post-pregnancy
- · Phlebitis, arteriosclerosis, ischemic vascular diseases, leg edema, pulmonary, edema.
- General immobility or bedridden situations post-operation; venous insufficiency
- Phoenix embolism stockings accelerate venous return flow rate and avoid blood clot formation which may cause pulmonary embolism or chronic venous deficiency
- Made of nylon and elastin material
- Contains no latex

ISO 13485 ISO 9001

[ | ISO 13485 | ISO 9001

ISO 13485 | ISO 9001

ISO 13485 | ISO 9001

### MEDI PHOENIX



### VARICOSE SOCKS (SOOTHING) 18-22 mmHg

Medi –Phonex soothing stockings are very thin, transparent, pain knit, soft and comfortable

Assists in graded pressure blood circulation and restricts problems caused by vessel failure

- Fatigueless and pain on legs, early venous deficiencies, or varicose stages where no edema is present yet, spider vans, pregnancy varicose vans, sclerosis of varicose vans
- CCLI 18/22 mmhg is the pressure class equivalent up to 24/28 kPa
- Medi –Phonex low pressure varicose stockings assist in avoiding such discomforts

#### Available in ;

- Knee High Varicose Stocking (Closed Toe)
- Thigh High Varicose Stocking (Closed Toe)
- Pantyhose Varicose Stocking (Closed Toe)
- Pregnancy Pantyhose Varicose Stocking (Closed Toe)

### MEDI PHOENIX



### VARICOSE SOCKS (SOOTHING) 20-30 mmHg

- Recommended for treatment of vascular diseases, severe symptoms, significant varicose vans or intermediate edema, upon elution of minor ulcers, post-operation for sclerosis and varicose vans, deep vein thrombosis and deep van deficiencies
- CCLII 20/30 mmHg is the pressure class equivalent up to 32/43 kPa.
- For diabetics with bunion and malnutrition of skin, "open Toe" option is available.

#### Available in:

- Knee High Varicose Stocking (Closed & Open Toe)
- Thigh High Varicose Stocking (Closed & Open Toe)
- Pantyhose Varicose Stocking (Closed & Open Toe)

# MEDI



### VARICOSE SOCKS (SOOTHING) 30-40 mmHg

- Constitution arising from constitutional or post-thrombotic venous insufficiency
- Strong tendency to edema
- Secondary varicose
- White atrophy, dermatosclerosis
- Following treatment for severe ulcers, if they have caused a relapse

#### Available in:

- Knee High Varicose Stocking (Closed & Open Toe)
- Thigh High Varicose Stocking (Closed & Open Toe)
- Pantyhose Varicose Stocking (Closed & Open Toe)

[ E | ISO 13485 | ISO 9001

[ [ 150 13485 | 150 9001

[ E ISO 13485 | ISO 9001

# **ORTHOPEDIC EQUIPMENTS**






# **ORTHOPEDIC EQUIPMENTS**



# Wicromed.



# Wicromed.



# $\mathcal{W}$ icromed.



### TRAVEL CUSHION

- Visco elastic
- Avoids neck cricks during travel and sleep
- İdeal for a comfortable and smooth sleep
- · Recommended for extended travels
- Occludes neck space assist in avoiding cervical disc hernia and spinal discomfort
- Zipper, removable and washable fabrics
- Fits snug around neck for maximum support
- Soft seat, transforms any hard surface into a comfortable seat

### **VISCO CUSHION**

- Cushions manufactured by Visco elastic technology are directly sensitive to body temperature, pressure and ambient temperature
- Fully supports neck space, and avoid chronic pains of neck spine articulation
- During sleep, ensures stabile line for the spine
- Hygienic
- Non-allergenic
- Manufactured free of CFC gas which is hazardous for human health
- Contains no carcinogen

Available
size
S
M
L

ISO 13485 | ISO 9001

### **LUMBAR CUSHION**

- Durable fabric covering and polyurethane foam filled
- Cushion helps the lumbar and sacral region of the spinal column
- A great carry along for any seat at the office, in the car, or a airplane ride
- Ensures correct position sitting thanks to anatomical form thereof
- · Readily removable and washable cover
- Considerably diminishes waist- and back pains
- Directly supports waist and vertebra
- Fits every style of seat to the elastic band at back
- Hygienic
- Manufactured free of CFC gas which is hazardous for human health
- · Contains no carcinogen
- Types; Visco and Polyurethane






# Wicromed.



### **DONUT CUSHION**

- Elevate the back of the rear for maximum weight relief
- Cushion made of polyurethane foam which is continuously elastic and may maintain form thereof
- Reduces hemorrhoids decubitus scars, hip and bone pain, and coccyx degeneration pain
- Reduce pressure points in the legs, upper thighs, and rear.
- Hygienic
- Manufactured free of CFC gas which is hazardous for human health
- · Contains no carcinogen
- Offer balance so there's no tilting when rotating or reaching. (This type of tilting can cause additional pain.)
- Improve posture, which can increase blood flow and circulation to encourage natural body healing.
- · Allow for easy travel for donut comfort when driving, working, sitting at home, and everywhere in between.

# **RESPIRATORY CARE**





# **RESPIRATORY CARE**









EOLO is a piston compressor nebulizer for aerosoltherapy. Designed for home-care with convenient handle and a large compartment for easy transport of unit, accessories and medications. transport of unit, accessories and medications.

Compressor: Olless and maintenance-free piston pump

Power Feeding: 230 V-50 Hz

Power consumption: 170 VA

Max Pressure: 2,5 bar 36 psi 250 Kpa

Operating Pressure: 1,1 bar

- Operating Pressure: 1, 10al

  Max Air Flow: 14 L/min

  Operating Air Flow: 5 L/min

  Neb-rate with 4 ml solution: 0,40 ml/min

  Mass Median Aerodinamic Diameter: 2,44 µm (MMAD)

  Respirable Fraction (< 5µm): >80 %

- Duty Cycle: Continuo
   Noise Level (EN 13544-1 Standard): 55 dB
- Weight: 1,65 Kg
   Size (Main Unit): 14x12x22 cm

#### included accessories;

- HI-FLO Jet Nebulizer
- · Adult Masks with elastic strings
- Child Masks with elastic strings
- Mouthpiece
- Nosepiece (non invasive)
- Air Tube
- Air Filter
- Carrying Bag







### RINA VAP COMPRESSOR NEBULIZER

RINA VAP is a piston compressor nebulizer for aerosoltherapy. Indicated for home-care and intensive use, its design includes a convenient handle for easy transport of unit and a large compartment for storage of accessories and medications.

- Power: 230 V AC /50HZ or 120V AC /60HZ
- Current:0.70 or 1.5A
- Max capicity of medicine cup: 6ml
- Particle size: 0.5 to 10um
- Average Nebulizer rate: ≥0.33 ml/min
- MMAD:3 0~6 0µm
- Noise: ≤ 60dBA
- Unit size: 270 x 200 x 110 MM
- Unit weight: 2.0 KG
- Suggested applications: hospital, clinical, nursing home, homeca-
- Accessory:
- 1.5M air tube
- Mouthpiece
- Nosepiece
- Adult mask · Child mask
- Air filter (5 pcs)
- Medicine cup







### **AEROSOL CHAMBER**

- Made of liquid silicone.
- Non-allergenic.
- Fully face fitting.
- Accelerates lung function progress.
- Avoids medicine loss of time delay during inhalation.

Size
S
M
1







### **OXYGEN CONCENTRATOR 7F-5**

Oxygen density: ≥ 90%Oxygen flux: 1~5L/min

• Output pressure: 0.04~0.06Mpa

Overrunning pressure alarm

• Nebulization rate: 0.2 ml/ dk

• Noise:  $\leq 54dB(A)$ 

• **Power supply:** AC220± 22V, 50± 1Hz

Input power: ≤ 500VADimensions: 380x310x600

• Weight: 25 kg

# Wicromed.



### **INCENTIVE SPIROMETERS**

The patient should lift two or three balls by deep inspiration, following the instructions of the doctor specialist, holding the inspiration as long as possible and removing the mouthpiece from his mouth to exhale. Can be fully disassembled for cleaning operations in between patients.

- 600-900-1200 ml
- Corrugated tube
- Mouthpiece
- Air filter
- Three balls

# **PERSONAL CARE PRODUCTS**





# **PERSONAL CARE PRODUCTS**



### Oncomed



### **SCALE -HEIGHT METER DIGITAL - SC 105**

- Maximum weight: 200 kg • Minimum weight: 2 kg
- Weighing accuracy: 100 gr
- Measurement range: 70-190 cm • Measurement accuracy: 0.5 cm
- Platform dimensions: 372x272 mm
- General dimensions: 535x275x940 mm
- Fixed weight: 13 kg
- Operating ambient temperature: 5 °C.-40 °C
- Humidity: <85 RH
- Power adapter: AC. 220V. Battery included; chargeable
- Digital display

### Oncomed



### SCALE - HEIGHT METER MECHANICAL -SC 101

- Maximum Weight: 120 kg. (265 liber)
- Lower sections of measurement indicator: 0.5 kg
- Height measurement range: 700-1900 mm. (271/2-491/2)
- Lower sections of measurement: 5 mm
- Platform dimensions (width x length): 375x270
- Device dimensions (WxLxH): 695x286x930 mm
- Weight: 25 kg.
- Mechanical display

### Oncomed



### TEST HCG - PREGNANCY TEST CARD

- · Accurate and quick result
- Solely for in-vitro diagnosis
- Disposable
- Cut off value: in urine 20 mlU/ml hCG
- Sensitivity: % 100 • Specificity: % 99, 9
- CE
- Made in Turkey

### Package includes:

- HCG pregnancy test device,
- Disposable plastic dropper.
- User's manual
- · Silica gel sachet as dehumidifier

 $\epsilon$ 

 $\epsilon$ 



ISO 15197 | ISO 9001

# Wicromed.



# Wicromed



### **COLD &HOT PACK (THERMO GEL)**

- Cold/Hot use
- Special design for painful sites of body
- Resistant to microwave and fridge
- Ergonomic design
- Made in Turkey

### Size

13x13 cm 13x26 cm 26x37 cm

### COLD &HOT PACK (THERMO GEL)

- Cold/Hot use
- Special design for painful sites of body
- Resistant to microwave and fridge
- Ergonomic design
- With Case
- Made in Turkey

### Size

13x13 cm


13x26 cm 26x37 cm

# **DISPOSABLES**





# **DISPOSABLES**









### **ORTHOPEDIC CAST (AMERICAN)**

- Made of synthetic knit by polyester or fiberglass material which activates by water and immersed in polyurethane resin
- Once removed does not cause allergy, itching or dust
- High X-ray permeability
- Skin perfectly breathes thanks to micro vents, hygienic
- Solid, resistant, lightweight and excellent hardness
- Readily applicable by length of 4,0 meters instead of 3,6

Size	
2 Inc	
3 Inc	
4 Inc	
5 Inc	

CE

### I.V CANNULA

- FEP (flour ethylene propylene) (Teflon) coated
- Tap port
- Low pressure valve system at port side
- · Readily open able injection cover
- Angular and grooved flaps
- Lure-Lock adjustmentTransparent flash back
- Radio opaque feature against X-rays
- Dia. and color code in accordance with worldwide standards
- Non-pyrogenic and non-toxic
- CE. ISO certified

Height (cm)	Colour
14G	Orange
16G	Grey
18G	Green
20G	Pink
22G	Blue
24G	Yellow
26G	Purple

[ ISO 13485 | ISO 9001

### **ENEMA KIT**

- From soft PVC materials
- 750 ml capacity
- Ease of use rectum tip
- Soft PVC material
- 110 cm kink resistance tubing
- Castile Soap
- Drape
- PE gloves
- Lubricant gel



# Wicromed.



# Wicromed.



# Wicromed.



# Wicromed.



### **GAUZE HYDROPHILE**

- 20 fibers white, hydrophile
- 100% cotton yarn
- 30x80 cm
- Sterile independent packaging
- Does not contain latex.

### **GAUZE COMPRESS**

- 20 fibers white, hydrophile
- 7.5x7.5 cm. 4 layers folding
- Made of 100% cotton yarn
- Gamma sterilization
- Packaging: packages of 25, 50 and 100

### **ELASTIC BANDAGE**

- Made of cotton and latex yarn
- Air-permeable texture
- 2, 5 fold extension once stretched
- · Skin-compliant color
- Cellophane packaging
- W/wo Packaging

Size
6cm x 1,5m
8cm x 1,5m
10cm x 1,5m
12cm x 1,5m
15cm x 1,5m
20cm x 1,5m

### **GAUZE BANDAGE**

- 20 fibers hydrophilic
- Made of 100% cotton yarn
- High absorbent against bleeding by virtue of hydrophilic feature
- Edges knitted
- Sterilizable

S	ľ	ze –
5cm	Χ	1001
7cm	Χ	100r
10cm	Х	100

15cm x 100m 20cm x 100m

[ E ISO 13485 | ISO 9001

**( ( I I SO 13485 | I SO 9001** 

 $\epsilon$ 

 $\epsilon$ 

# **DISPOSABLES**

# **Wicro**Plast



### **Wicro**Pore



### **Wicro**Trans



### **WicroFix**



### **COTTON FABRIC PLASTER**

- Zinc oxide rubber based adhesive medical fabric.
- · May be torn off properly in length and breadth.
- No need to use scissors.
- Used for fixing purposes in all kinds of surgical application.
- Remains adhesive for a long period.
- · Not visible through X-ray photoraph.
- Made in Turkey

#### Available size

5m x 1.25cm 5m x 2.5cm 5m x 5cm 5m x 7.5cm 5m x 10cm

**( (** 150 13485 | 150 9001

### NON-WOVEN HYPOALLERGENIC PLASTER

- Transtsmits air and sweat but water due to its micro pores. Its effects to the skin, even to very sensitive skin are very little due to its special adhesive. Sticks easily on hairy areas, does not cause pain when pulling out.
- No need to use scissors becuse of its special perforated area.
- May be torn off properly in length and breadth.
- Not visible through X-ray photoraph.
- Made in Turkey

### POLYETHYLENE HYPOALLERGENIC PLASTER

- Trantsmits air and sweat due to its micro pores, water resists
- Its effects to the skin, even to very sensitive skin are very little due to its special adhesive.
- Can be torn off from length and breadth easily.
- Remains adhesive for a long period.
- Application area may be seen due to its transparent surface.
- Not visible through X-ray photoraph.
- Made in Turkey

### HYPOALLERGENIC ELASTIC PLASTER

- · Well-adjusting to the body due to its elastic feature.
- Convenience on areas with edema and on moving areas(orthopedic,abdolminal,thorax obstetrics)on all kinds of surgical attempt areas, on shoulders, elbows, armpits, process areas, in post-operative usage.
- Can be pulled out cleanly and painlessly.
- Made in Turkey







ISO 13485 ISO 9001

# **Wicrolpek**



# Wicromed.



# **Wicromed**.







### **SILK PLASTER**

- Acrylic on Zinc oxide,rubber based,adhesive,medical silken plaster made of acetate silk.
- May be torn off propely in length and breadth. No need to use scissors.
- Used for fixing purposes in all kinds of surgical application.
- Remains adhesive for a long period.
- Not visible through X-ray photograph.
- Made in Turkey

### Available size

5m x 1.25 cm 5m x 2.5 cm 5m x 5 cm 5m x 7.5 cm

5m x 10cm

ISO 13485 | ISO 9001

### **ADHESIVE EYE PAD (ELIPTIC)** (FOR CHILDREN)

- Hypoallergenic –adhesiveness suitables for skin-sense-
- Light-proof black pad provides the dark environment needed for eye
- Designed for children; safety and comfortable use after eye injuries and postoperative eye cares
- Steril until the pack is opened
- Can be used on both two eyes and behind glasses
- Single use only
- Made in Turkey

### STERILE URINE BAG

- 2000 ml. volume
- Made of medical PVC
- Hook hole
- Tube made of PVC material of 90 cm
- Appropriate hardness which avoids bending
- W/wo tap
- Disposable and sterile

### PEDIATRIC STERILE URINE BAG

- To collect urine sample from infants
- 10 ml level indicators on bag
- Made of medical PVC
- Independent and sterile
- Girl and boy option







 $\epsilon$ 

# **DISPOSABLES**

Aloe Mate



Oncomed



Wicromed.



Wicromed.



### **EXAMINATION GLOVES WITH ALOE VERA**

- Inner surface of gloves is coated by natural, pure, superior quality and soft aloe vera gel.
- Activates by body temprature and ensures moisture of hands by gel, avoids dehydration.
- Contains no latex proteins.
- Powder free
- Single use only

Packaging:

Box:100 units Package: 20 boxes

Size
S
M
L
XL

**SCALP VEIN SET** 

- Cannula: Thin walled stainless needle
- Sterile Non toxic Pyrogenic free
- Wing: Color coded flexible wing for easy manipulation
- Tube Capacity: 0.43mL for 300mm tube length
- Connector: Lock hard connector (luer lock)
- Soft connector (luer slip) available upon request
- Tube Length: 12" (30 cm), available upon request

Gauge (mm)	Needle O.D
19	1.0
20	0.9
21	0.8
23	0.6
25	0.5

**BOUFFANT CAP** 

- Caps comply with NFPA standards and are latex-free.
- Double-stitched elastic band keeps caps flat and comfortable on the forehead.
- · Economy cap, lint free.
- Made of spunbond polypropylene.
- Antiallergenic
- · Colour: White and Green
- Packaging:
   Box:100 units
   Package: 250 boxes
- Made in Turkey

**ELASTIC MASK - 2 LAYERS** 

- Hygienic and anti-allergenic
- Air permeable ,liquid -tight
- Ultrasonic suture
- Head and beard covers provide complete coverage for worry-free performance during surgery
- Available in a spunbound polypropylene or multi-laver material
- Disposable
- Single use only
- Packaging: 150 units
- Made in Turkey

 $\epsilon$ 

ISO 13485 | ISO 9001

 $\epsilon$ 

 $\epsilon$ 



### **ELASTIC MASK - 3 LAYERS**

- · Hygienic and anti-allergenic
- Air permeable ,liquid -tight
- Ultrasonic suture
- Head and beard covers provide complete coverage for worry-free performance during surgery
- Available in a spunbound polypropylene or multi-layer material
- Disposable
- Single use only
- Packaing: 50 units
- Made in Turkey



### **NASAL OXYGEN CANNULA ADULT / PEDIATRIC**

- Made of non-toxic ,non-irritating PVC
- 2 m length.
- Soft and clear
- Air passage saffety
- Soft cannula resistant to bending
- Clamp for narrowing and expansion
- Single use only
- Made in Turkey



### **OXYGEN MASK - ADULT / PEDIATRIC**

- Made of medical PVC material
- 2 m.hose
- Adjustable nose clip assures comfortable fit
- Standard connection connector
- Elatic rubber (ended)
- Air passage safety
- · Clear, soft vinyl for patient comfort and visual assessment
- Made in Turkey



### **OXYGEN MASK (with CHAMBER) ADULT / PEDIATRIC**

- Made of medical PVC material
- 8 ml. Medicine chamber
- 2 m.hose
- · Adjustable nose clip assures comfortable fit
- Standard connection connector
- Elatic rubber (ended)
- Air passage safety
- Clear, soft vinyl for patient comfort and visual assessment
- Made in Turkey

[ ISO 13485 | ISO 9001


ISO 13485 | ISO 9001

[ [50 13485 | 150 9001



[ ISO 13485 | ISO 9001

# **DISPOSABLES**



ISO 13485 | ISO 9001







# Wicromed.



### **FOLEY CATHETER SILICONE**

- Made of Non-toxic, non-irritating silicone material
- X-ray line
- · Two lateral aperture on distal end
- Two-ways
- Vacuum control
- Disposable and individually packaged
- **Size:** 6FR 24FR

[ ISO 13485 | ISO 9001

### **EXTERNAL MALE CONDOM CATHETER**

- Made of natural latex
- Proximal end designed to ensure easy connection with urine bag
- Self-adhering band
- Disposable and individually packaged

### **ULTRASOUND - ECG GEL**

- · Water-soluble, easily cleaned
- Hypoallergenic
- Free from fatty substance, formaldehyde and salt
- Non-toxic and odorless
- May be used in Ultrasound and ECG
- PH 7
- Manufacture under vacuum
- Free from air bubble
- Transparent in blue color
- 1000 ml
- Made in Turkey













### **FOR A**



### **ORAL CARE SET**

- Single packed, EO Sterile
- 100 ml. solution
- 20 Medical Sponge swaps easily cleaned
- Moisturizing Gel and scaled container
- Mouth Antiseptic
- Hydrogen peroxide %1, 5

### HYDROPHILIC CATHETER WITH **WATER SOCKET**

- Maintains slickness for extended periods once get into contact with water
- Avoids trauma by virtue of rounded end
- Unrecognizable precise drainage holes ensures easy discharge of bladder
- Non-traumatic non-infection application thanks to nontoxic gel on surface
- · Low-friction catheter minimizes risk of infection, while maximizing comfort
- Sterile

MALE	FEMALE
06	06
80	08
10	10
12	12
14	14
16	16
18	

ISO 13485 | ISO 9001

### **DIGITAL THERMOMETER (FLEX)**

- Digital LCD screen
- Measuring method: In the mouth (oral), In the anus ( rectal )
- · Compact And Sensitity
- Range: 32.0-42.0 °C
- Automatic Shut Down
- · Very Light and user friendly
- Size:12mm x 22mm x 145mm
- Storage of the last measured value
- 100% waterproof
- · Without glass, mercury free
- · Automatic switch off
- Large Display
- · Can be disinfected
- · High accuracy

ISO 13485 | ISO 9001

#### DIGITAL THERMOMETER IR16 EAR FORFHFAD

- **Dimensions:** 120mm (L) x 15mm (W) x 25mm (H)
- Weight: 69g
- Precision: 0.1°C
- Accuracy: 1.36.0°C 39.0°C range ±0.2°C
   2 <36.0°C and <39.0°C range ±0.3°C

   Power supply: 1 units of CR2032 3V lithium battery
   Battery life: At least 1000 measurements

- Memory: Latest 10 measurements
   Power protection: 30 seconds
- Storage conditions: 5°C 40°C range, 95% or less relative humidity
- Transportation conditions: 20°C 50°C range 95% or less humidity
- Measurement range: 0°C 50°C
   Measurement unit: °C or °F
- Measurement site: Ear and Forehead mode
- High/Low measurement:


0°C LCD Display LO 0°C~37.9°C LCD Normal 38.0°C~50°C LCD Normal / HI 50°C LCD Display HI

- Lens Filter/Magnifier: Yes
- Self-Diagnosis Function: Warning Signal of Abnormal Conditions to Prevent Inaccurate Measurement



ISO 13485 | ISO 9001





# <mark>mt</mark> sağlık ürünleri